



Creating a Cancer-free World. One Person, One Discovery at a Time.

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

## FuSpot: A Web-based Tool for Visual Evaluation of Fusion Candidates

**Jackson A. Killian**  
Oct. 24, 2017

# Outline

- Fusions and Cancer
- Fusion Detection
- FuSpot
- Methodology
- Case Study


The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

# Gene fusions are intimately linked with cancer.


- First genetic defect linked with cancer was t(9;22)
  - Philadelphia Chromosome → CML
- Created hybrid gene BCR-ABL
  - Dysregulated tyrosine kinase producer
  - Rapid cell division, oncogenic



The James

# Studying fusions leads to better therapies.

- Chemotherapy ineffective, transplants dangerous
- Developed BCR-ABL specific drug
  - Tyrosine kinase inhibitor – Imatinib
- Long term survival increased from 30% to ~90%<sup>1</sup>
  - 98% complete response, 17% relapse rate
- Also discovering prognostic markers
  - *PAX-FKHR* risk group stratification in ARMS<sup>2</sup>
  - Similar implications for *TMPRSS2-ERG* in Prostate Carcinoma<sup>3</sup>



1. Druker BJ, et al. New England Journal of Medicine 355, 2408–2417 (2006)
2. Sorensen PH, et al. J Clin Oncol. 2002 Jun 1;20(11):2672-9.
3. Berg KD. Dan Med J. 2016 Dec;63(12). pii: B5319.

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

# Many fusion detectors are available for download.

- Therapy implications + advent of RNA-seq led to development of many sequencing-based fusion detection tools ( $\geq 30$  unique)<sup>4,5</sup>
- Some emphasize specificity, but sacrifice sensitivity
  - Trusted results, but miss true positives
- Others trade specificity to ensure good sensitivity
  - Large numbers of False Positives
  - Time consuming and costly PCR validation

4. Liu et al. Nucleic Acids Res. 2016; doi:10.1093/nar/gkv1234.

5. Kumar et al. Wiley Interdiscip Rev RNA. 2016 Nov;7(6):811-823. doi: 10.1002/wrna.1382

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

# Our tool, FuSpot, accelerates discovery.


- **The Problem:** Difficult to choose tool that will get all true positives without overload of false candidates
- **Solution:** Use a sensitive tool + our tool **FuSpot** to filter out false positives with little time and no experimental cost

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

# FuSpot has an intuitive workflow.



The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

## Breakpoint Coordinates:

References: ☒ Auto ☐ Custom

Enter the genomic coordinates of a fusion breakpoint and FuSpot will retrieve the genomic sequences as well as the sequences of the nearest exons of both fusion gene partners to use for alignment:

Genome Build:

### 5' End

Coordinate:

Strand:

Reference Length:

Gene Name:

### 3' End

Coordinate:

Strand:

Reference Length:

Gene Name:

Title:

Alignment Type:

[Generate Example Breakpoint](#)

\*All file sizes must be less than 4Mb.

## Single-End Reads:

[Read File\\*](#)

## Paired-End Reads:

[First Mate File\\*](#)

[Second Mate File\\*](#)

[Get Example Fasta](#)

## Align and Visualize

[Generate](#)

# FuSpot Interface

## Inputs

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER



## Breakpoint Coordinates:

References: ☒ Auto ☐ Custom

Enter the genomic coordinates of a fusion breakpoint and FuSpot will retrieve the genomic sequences as well as the sequences of the nearest exons of both fusion gene partners to use for alignment:

Genome Build:

### 5' End

Coordinate:

Strand:

Reference Length:

Gene Name:

### 3' End

Coordinate:

Strand:

Reference Length:

Gene Name:

Title:

Alignment Type:

[Generate Example Breakpoint](#)

\*All file sizes must be less than 4Mb.

## Single-End Reads:

[Read File\\*](#)

## Paired-End Reads:

[First Mate File\\*](#)

[Second Mate File\\*](#)

[Get Example Fasta](#)

## Align and Visualize

[Generate](#)

# FuSpot Interface

## Inputs

1. Genomic Coordinates of Breakpoint

Ex.

Chr2:35479453 +

Chr17:27374426 -

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

## Breakpoint Coordinates:

References: ☒ Auto ☐ Custom

Enter the genomic coordinates of a fusion breakpoint and FuSpot will retrieve the genomic sequences as well as the sequences of the nearest exons of both fusion gene partners to use for alignment:

Genome Build:

### 5' End

Coordinate:

Strand:

Reference Length:

Gene Name:

### 3' End

Coordinate:

Strand:

Reference Length:

Gene Name:

Title:

Alignment Type:

[Generate Example Breakpoint](#)

\*All file sizes must be less than 4Mb.

### Single-End Reads:

[Read File\\*](#)

### Paired-End Reads:

[First Mate File\\*](#)

[Second Mate File\\*](#)

[Get Example Fasta](#)

### Align and Visualize

[Generate](#)

# FuSpot Interface

## Inputs

1. Genomic Coordinates of Breakpoint
2. Reads adjacent to the breakpoint\*

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

Strand: +/-

Reference Length: 10-1000

Gene Name: Optional Gene Name

Alignment Type: Local

Generate Example Breakpoint

\*All file sizes must be less than 4Mb.

**Single-End Reads:**

Read File\*

**Paired-End Reads:**

First Mate File\* Second Mate File\*

Get Example Fasta

**Align and Visualize**

Generate

# FuSpot Interface

## Inputs

1. Genomic Coordinates of Breakpoint
2. Reads adjacent to the breakpoint\*

\*FuSpot provides extraction tool



### How to get reads to input to FuSpot?

Align with a Chimeric-capable RNA aligner. Then download our read extraction tool and run it on your aligned file to get reads local to the candidate fusion.

Download Extraction Tool

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

## Breakpoint Coordinates:

References: ☒ Auto ☐ Custom

Enter the genomic coordinates of a fusion breakpoint and FuSpot will retrieve the genomic sequences as well as the sequences of the nearest exons of both fusion gene partners to use for alignment:

Genome Build:

### 5' End

Coordinate:

Strand:

Reference Length:

Gene Name:

### 3' End

Coordinate:

Strand:

Reference Length:

Gene Name:

Title:

Alignment Type:

[Generate Example Breakpoint](#)

\*All file sizes must be less than 4Mb.

### Single-End Reads:

[Read File\\*](#)

### Paired-End Reads:

[First Mate File\\*](#)

[Second Mate File\\*](#)

[Get Example Fasta](#)

### Align and Visualize

[Generate](#)

# FuSpot Interface

## Inputs

1. Genomic Coordinates of Breakpoint
2. Reads adjacent to the breakpoint\*

\*FuSpot provides extraction tool

## Align and Visualize


The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER


# Reference Extraction

- Gather sequences of each gene partner
- Collect bases outward in both directions from breakpoint
  - Genomic
  - Exonic
- UCSC canonical exons



# FuSpot Alignment

- References:
  - Automatic Retrieval
    - Human or Mouse
  - Custom
    - Any organism
    - Any combination of gene transcripts
- Simultaneous alignment to each reference
- Determine right-to-left connection (if any)



The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

# Case Study: Edgren et al.<sup>6</sup> and BEERS<sup>7</sup>

- Edgren et al. contains several well known fusions
  - True Positive case
- BEERS is synthetic dataset containing no fusions
  - False Positive case
- Ran FusionCatcher and FusionMap on both sets


6. Edgren et al. Genome Biology. 2011;12(1). doi:10.1186/gb-2011-12-1-r6.

7. Grant et al. Bioinformatics. 2011. doi:10.1093/bioinformatics/btr427.

The James




THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER



## True Positive Case


- Reported by both tools
- Support breakpoint
  - Group 1+4
- Non-fusion Reads
  - Group 2+3
- Well supported breakpoint




The James



# False Positive Case




- Reported by FusionMap
  - Claims 9 supporting reads
- Non-fusion reads
  - Group 1
- Seem to support breakpoint
  - Group 2
- No biological sense
  - Group 3
  - Suggests homology between fusion partners



The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER



## Realign Groups 2+3

- Align using only 5' partner references
- All reads align → Sequence homology!
- This is what tricked the detector
- Confidently categorize as False Positive

The James

# FuSpot solves the tradeoff problem.

- Fusion Detectors vary in tradeoff between specificity/sensitivity
  - Risk missing important fusions
  - Or spend resources validating false positives
- FuSpot visualizes evidence in easily digestible form
- Enables convenient and rapid inspection by researcher

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER

# FuSpot is online and facilitating research.

- FuSpot takes advantage of human reasoning to eliminate tedious post-processing
- FuSpot will result in more fusions to be identified at lower overall cost
- Already supporting research in The James
  - Dr. He and Dr. de la Chapelle
  - Dr. Cynthia Timmers, Director of The Solid Tumor Translational Science Shared Resource at OSUCCC  
The James

# Thank You

PELOTONIA® 17

Authors:

Jackson A Killian, Taha M Topiwala, Alex R Pelletier,  
David E Frankhouser, Pearly S Yan, Ralf Bundschuh

<http://bioserv.mps.ohio-state.edu/FuSpot/>

The James



THE OHIO STATE UNIVERSITY  
WEXNER MEDICAL CENTER